

Infant Radiant Warmer



STANDARD COMPOSITION:

- Temperature Control Mode: pre-heating control, servo control & manual control;
- Instant Parameter Display: set temp, skin temp, heat power rate;
- Built-in APGAR Timer and $>37^{\circ}\text{C}$ preset function
- Audible/Visual Alarms: temp deviation, sensor failure, power failure, temp override;
- Tilting Bed and Rotary Radiant Head;
- X-ray Cassette Tray under infant bed, examination light of high illuminance;
- Use anti-blast quartz infrared tube as radiant heating source.



CCI-TD1001

Phototherapy Lights
Mounting beside radiant head
of Infant Radiant Warmer.

Optional Devices (At Buyers' Selection):



CCI-1003

Intra Low-pressure suction.



DG

Storage drawer



SJ

Electric height adjustable

Specifications:

Specifications	CCI-1000	CCI-2000	CCI-2000(B)
Operational condition:			
Environment temperature	+18°C ~ +30°C	+18°C ~ +30°C	+18°C ~ +30°C
Environment relative humidity	30% ~ 75%	30% ~ 75%	30% ~ 75%
Environment air velocity of flow	<0.3m/s	<0.3m/s	<0.3m/s
Power requirement	≤ 1000VA	≤ 1000VA	≤ 1000VA
Temperature control mode	Auto	Pre-heat, auto and manual	Pre-heat, auto and manual
Temperature Control Range	25°C ~ 37°C	25°C ~ 37°C	25°C ~ 37°C
Temperature override mode to	-----	38°C	38°C
Temperature sensor accuracy	≤ 0.3°C	≤ 0.3°C	≤ 0.3°C
Mattress temperature uniformity	≤ 2.0°C	≤ 2.0°C	≤ 2.0°C
Temperature rise time	≤ 45 min	≤ 45 min	≤ 45 min
Heat output indicate	0~100%	0~100%	0~100%
Radiant box horizon angle adjust range	0~±90° step-less	0~±90° step-less	0~±90° step-less
Mattress Tray Size	68cm×56cm	68cm×56cm	68cm×56cm
Mattress tilt adjust range	0~±10° step-less	0~±10° step-less	0~±10° step-less
Mattress tray Max. load	10 kg	10 kg	10 kg
Transfusion pole Max. load	2 kg	2 kg	2 kg
Shelf Max. load	2 kg	2 kg	2 kg
X-ray	Yes	Yes	Yes
RS-232 Connector	----	Yes	Yes
Mattress to Radiant box	81 cm	81 cm	81 cm
Heater infrared wavelength	1 ~ 3 μm	1 ~ 3 μm	1 ~ 3 μm
APGAR Timer	-----	Yes	Yes
Alarms			
Power failure alarm	Yes	Yes	Yes
Sensor failure alarm	Yes	Yes	Yes
Temperature deviation alarm	±1.0°C	±1.0°C	±1.0°C
Over temperature alarm	≤ 38.0°C	≤ 38.0°C、≤ 39.0°C	≤ 38.0°C、≤ 39.0°C
System failure alarm	-----	Yes	Yes
Low-pressure suction	-----	-----	Yes
Pressure adjust range	-----	----	0 ~ 22Kp(0 ~ 165mmHg)
Operating noise	-----	----	≤ 55dB
White Phototherapy light	-----	-----	1 set
7L Oxygen cylinder (with oxygen manometer & oxygen flow meter)	-----	-----	2 sets
Neonate resuscitation bag	-----	-----	1 set
Packing size	110×90×110 cm	110 ×90 ×110 cm	110 ×90 ×110 cm
Net weight	85 kg	85 kg	95 kg
Gross Weight	115 kg	115 kg	125 kg

China Care Medical Equipment Co., Ltd

Address: DaXin Building 2203C room, DeZheng North Road No. 538, YueXiu District, Guangzhou, China. 510045

Cellphone/Whatsapp: 008613926480041 **Skype:** chinacaremed **Wechat:** Joseph20130708

Facebook/Youtube/Instragram/LinkedIn: china care medical

Website: www.chinacaremedical.com **Email:** chinacaremed@gmail.com