

Product Introduction

Model : CCE-3303B

Product Description



Product Features :

1. Portable and delicate design, easy to operate
2. Accurate pulse pace identification function
3. High accuracy digital filter, automatic baseline adjustment
4. Working modes: Manual, Automatic, Arrhythmia Analysis, Storage
5. 80mm, three channel format recording, excellent automatic interpretation
6. 320x240 graphic 3.5 inch color LCD to simultaneously display ECG information
7. 250 patient cases storage and replay (SD card is optional)
8. Detailed patient information record
9. Freeze function
10. Adapt to 110-230V,50/60Hz power supply. Built-in rechargeable Ni-MH battery
11. RS232 and USB interface, optional for USB storage, external laser printer printing, PC ECG software functions.

Technical Specification :

Lead	Standard 12 leads	Time constant	$\geq 3.2s$
Lead acquisition	12bit/1000Hz (synchronously 12 leads)	Frequency response	0.05~160Hz(-3db)
Work mode	Man / Auto/Analysis/Storage	Noise level	$< 15\mu Vp-p$
Filter	AC Filter : 50Hz / 60Hz EMG Filter : 25Hz / 45Hz Anti-Drift Filter : 0.15Hz (Adaptive)	Inter-channel Interference	$\leq 0.5mm$
CMMR	$\geq 120dB$ (with AC FILTER)	Sensitivity	5, 10, 20, 40, Auto (mm/mV)
Input circuit	Floating; Protection circuit against Defibrillator effect	Recording mode	1ch+, 3ch, 3ch+ Default Man. Mode is 3ch format Default Auto. Mode is 3ch+ format
Input Impedance	$\geq 50M\Omega$	Paper speed	6.25、12.5、25、50mm/s ($\pm 3\%$)
Input circuit current	$\leq 0.05\mu A$	Paper	80mm *20m high resolution thermal roll paper
Patient current leakage	$< 10\mu A$	Communication	RS232/ USB port
Calibrating voltage	1mV $\pm 2\%$	Power Supply	AC 110-230V($\pm 10\%$),50/60Hz ($\pm 1Hz$),60VA DC Built-in rechargeable Ni-MH battery,12V (1500mAh)
Voltage tolerance	$\geq \pm 500mV$	Dimension/Weight	300mm \times 230mm \times 75mm / 2.8Kg

Packing List :

1. Main Unit : 1 set
2. Patient Cable : 1pc
3. Chest Electrodes : 1set (6 pcs)
4. Limb Electrodes : 1set (4 pcs)
5. Power Cable : 1pc
6. Grounding Cable : 1pc
7. Fuse : 2pcs
8. Paper Axis : 3pcs
9. Thermal Paper : 1 roll
10. Use Manual : 1 copy
11. The factory Inspection Report : 1 copy

Picture for Accessories :

